Opioid-Substitutions-
therapeutika in Vollblut

Werte-Ermittlung: OSD RV 25-10

Opioid Substitution Drugs
in Whole Blood

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1) The target value is calculated as the average of multi-
ple analyses by independent laboratories accredited
according to DIN EN ISO 17025 resp. DIN EN ISO
15189.
The analyses were performed by GC/MS and
(UP)LC/MS(-MS) according to the guidelines of the
German Medical Association for Quality Assurance
(RiliBÄK) and the guideline of the GTFCh (Society of
Toxicological and Forensic Chemistry) for quality as-
surance of forensic toxicological research [1] [2].

2) The estimated expanded uncertainty U of the target
value is calculated with a coverage factor of k=2 cor-
responding to a level of confidence of about 95% as
defined in the GUM [3].

3) Confidence Interval:
Range in which the real value is located with a signifi-
cance level of 95% (target value ± U, k=2).

4) The control ranges correspond to the accepted ranges
used in the proficiency tests of the GTFCh. They are
calculated from the target value and the standard de-
viation according to Horwitz (SD_{Horwitz}) [4] (last up-
date: 2024-06).
a) target value ± 30%
a) target value ± SD_{Horwitz}
Regardless, the control ranges should be defined by
the laboratory according to its own analytical require-
ments.

5) Density (20°C): 1.0523 ± 0.0001 g/cm³

6) The variation of the filling volume is max. ±0.5%.

7) This given value is only for your orientation. The data
required to establish a statistical evaluation were not
obtained.

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3) Guide to the Expression of Uncertainty in Measurement, ISO, 1995

4) HORWITZ, W. (1982): Evaluation of analytical methods used for regulation of foods and drugs. Anal Chem 54 (1): 67A-76A

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In-vitro-Diagnostikum

Opioid Substitution Drugs in Whole Blood

Description

Medidrug® OSD VB is a lyophilized reference material, prepared from human whole blood.

Intended Use

The reference material is used for quality assurance of the determination of opioids / opioid substitution drugs in whole blood.

Transport and Storage

The reference material can be shipped at ambient temperature. Upon receipt the material has to be stored at +2 to +8°C.

Stability

Originally closed: +2 to +8°C until the expiration date specified in the certificate of analysis.

After reconstitution: +2 to +8°C, 7 days

EDDP and 6-MAM have to be analysed on the day of reconstitution!

Preparation

Reconstitute each vial exactly with 2.0 ml of HPLC-grade or bidistilled water. Replace the stopper, swirl gently and allow the reference material to stand for 30 minutes at room temperature, then gently invert the vial for 15 minutes (rotation mixer) to ensure homogeneity. After reconstitution the reference material should be treated the same as patient specimens and run in accordance with the instructions accompanying the instrument, kit or reagent being used. Before sampling, allow the reference material to sit at room temperature and gently invert the vial for 10 minutes (rotation mixer) to ensure homogeneity.

Precautions

Each serum or blood donor unit used to manufacture this control was tested by state accepted testing methods and found non-reactive for following Antigens and Antibodies: HBsAG, anti-HIV-1, anti-HIV-2, anti-HBc and anti-HCV. No test method can guarantee complete assurance that products containing human source materials will be absent of these or other infectious agents. Therefore, handle them with the same precautions used with patient specimens.

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