Leervollblut (human)

Werte-Ermittlung: RV 20-18

Blank Whole Blood (human)

Value Determination: RV 20-18

Sang total témoin (humain)

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Notes

- 1) The target values were determined by independent accredited laboratories (EN ISO 17025 resp. DIN EN ISO 15189) according to the guidelines of the GTFCh (Society of Toxicological and Forensic Chemistry) and the German Medical Association for Quality Assurance (RiliBÄK). The target values were calculated by the median values of the Limit of Detection (LOD) of the respectively analytes [1][2]. Due to different capacity of the analytical methods by the laboratories, measured values may lie significantly beneath the denoted target values.
- 2) Method: A = UPLC-MS/MS; B = LC-MS/MS
C = GC-MS/MS; D = GC-MS
- 3) pH: 7.00
- 4) Density (20°C): 1.0521 ± 0.0001 g/cm³
- 5) Variation of the filling volume: max. ±0.5%.

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Blank whole blood (human)

Description

Medidrug Basis-line VB is a lyophilized reference material, prepared from human whole blood.

Intended use

For base-line determination in different applications or as negative control.

Check of the matrices

The matrices are analysed by independent institutions of forensic medicine and analytical laboratories.

The values listed apply only to the particular batch (Lot No. see label).

Transport and Storage

The reference material can be shipped at ambient temperature. Upon receipt the material has to be stored at +2 to +8°C.

Stability

Originally closed: +2 to +8°C, 36 months

But not longer than the expiration date specified in the certificate of analysis.

After reconstitution: +2 to +8°C, 30 days, stored tightly capped and dark

Preparation for use

Reconstitute each vial exactly with 2.5 ml of bidistilled or HPLC-water. Replace the stopper, swirl gently and allow the reference material to stand for 30 minutes at room temperature, then gently invert the vial for 15 minutes (rotation mixer) to ensure homogeneity.

After reconstitution the reference material should be treated the same as patient specimens and run in accordance with the instructions accompanying the instrument, kit or reagent being used. Before sampling, allow the reference material to sit at room temperature and gently invert the vial for 10 minutes (rotation mixer) to ensure homogeneity.

Precautions

Each serum or blood donor unit used to manufacture this control was tested by state-accepted testing methods and found non-reactive for following Antigens and Antibodies: HBsAG, anti-HIV-1, anti-HIV-2, anti-HBc and anti-HCV. No test method can offer complete assurance that products containing human source materials will be absent of these or other infectious agents. Therefore, handle them with the same precautions used with patient specimens.

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MEDICHEM Diagnostica GmbH