Medidrug® Ethanol S 500

8°C

2°C

LOT

CI702

2027-09

Edition 9 / 2026-02

Ethanol
in Serum (human)

Werte-Ermittlung: RV EtOH 13-05

Ethanol
in serum (human)

Determination of values: RV EtOH 13-05

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Détermination des valeurs: RV EtOH 13-05

Analysenmethode Method of Analysis Méthode d'analyse	Einheit Unit Unité	Verifizierter Wert Verified value Valeur vérifiée	Konfidenzbereiche / Confidence interval / Marge de confiance		
			I	II	III
GC Gaschromatographie / Gas chromatography / Chromatographie gazeuse	g/L	5,045	4,899 - 5,191	4,793 - 5,297	4,591 - 5,499
	g/kg	4,943	4,800 - 5,086	4,696 - 5,190	4,498 - 5,388
	mmol/L	109,51	106,34 - 112,68	104,03 - 114,98	99,65 - 119,36
ADH Alkohol-Dehydrogenase / Alcohol dehydrogen- ase / Alcool déhydrogénase	g/L	4,969	4,870 - 5,068	4,721 - 5,217	4,522 - 5,416
	g/kg	4,868	4,771 - 4,965	4,625 - 5,112	4,430 - 5,307
	mmol/L	107,86	105,71 - 110,01	102,46 - 113,25	98,15 - 117,56
GC + ADH Kombination / Combination	g/L	5,007	4,884 - 5,130	4,757 - 5,257	4,556 - 5,458
	g/kg	4,906	4,785 - 5,027	4,660 - 5,151	4,464 - 5,347
	mmol/L	108,68	106,01 - 111,35	103,25 - 114,12	98,90 - 118,46

Anmerkungen

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Notes

- 1) The verified value is calculated as the average of multiple analyses by independent laboratories accredited according to EN ISO 17025 resp. DIN EN ISO 15189.
The printed values correspond with serum.
- 2) Density (20°C): 1.0207 g/cm³
- 3) Confidence interval I – statistic:
Range in which the value is located with a significance level of 95%.
- 4) Confidence interval II – for forensic use^{a)}:
± 5% for ethanol concentrations ≥ 1.236 g/L
± 0.062 g/L for ethanol concentrations < 1.236 g/L
The measuring results must be rounded down after the second decimal place.
- 5) Confidence interval III – for clinical use^{b)}:
± 9% for ethanol concentrations 0.6 < [x] ≤ 5.0 g/L
± 15% for ethanol concentrations 0.2 ≤ [x] ≤ 0.6 g/L

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Ethanol in serum (human)

Description

Medidrug® Ethanol S is a liquid reference material, prepared from 100% human serum in 13 levels.

Intended use

Calibration or quality monitoring of determinations of ethanol in serum.

Assignment of Values

The verified value is calculated as the average of multiple analyses by independent accredited laboratories (EN ISO 17025 resp. DIN EN ISO 15189).

The indication of the confidence intervals for forensic and clinical applications are in accordance with relevant guidelines [1] [2].

The verified values listed apply only to the particular batch (LOT No. see label).

Transport and Storage

The reference material can be shipped at ambient temperature. Upon receipt the material has to be stored at +2 to +8°C.

Stability

Originally closed: +2 to +8°C, 36 months
But not longer than the expiration date specified in the certificate of analysis.

After working close the vials immediately and store them tightly capped at +2 to +8°C, do not let them standing open!

Ampoules are only for working at once!

Procedure

The reference material is ready to use and should be treated the same as patient specimens and run in accordance with the instructions accompanying the instrument.

Before sampling, allow the reference material to stand at room temperature and gently invert the vial to ensure homogeneity (5 min.).

Limitation of procedure

If there is evidence of microbial contamination or excessive turbidity in the used control, discard the vial.

Attention

Each serum or blood donor unit used to manufacture this control was tested by state-accepted testing methods and found non-reactive for the following Antigens and Antibodies: HBsAG, anti-HIV-1, anti-HIV-2, anti-HBc and anti-HCV. No test method can offer complete assurance that products containing human source materials will be absent of these or other infectious agents. Therefore, handle them with the same precautions used with patient specimens.

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MEDICHEM Diagnostica GmbH