Ethanol + Alkohole
Leervollblut (human)

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Ethanol + Alcohols
Blank whole blood (human)

Determination of target values: ALC 21-04

Ethanol + alcools
Sang total témoin (humain)

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Artikelnummer
Code
Code article

Ethanol + Alcohols
in Whole Blood (human)**Description**

Medidrug® ALC VB is a liquid reference material, prepared from 100% human blood in 10 levels.

Intended Use

Calibration or quality monitoring of determinations of ethanol, methanol, acetone and propan-2-ol in whole blood.

Assignment of Values

The values provided here in this insert were established by independent institutions working in forensic medicine.

The indication of the confidence intervals for forensic and clinical applications are in accordance with relevant guidelines (see bibliography).

The verified values listed apply only to the particular batch (LOT No. see label).

Transport and Storage

The reference material can be shipped at ambient temperature. Upon receipt the material has to be stored at +2 to +8°C.

Stability

Originally closed: +2 to +8°C, 36 months

But not longer than the expiration date specified in the certificate of analysis.

After working close the vials immediately and store them tightly capped at +2 to +8°C, do not let them standing open!

Ampoules are only for working at once!

Procedure

The reference material is ready to use and should be treated the same as patient specimens and run in accordance with the instructions accompanying the instrument.

Before sampling, allow the reference material to stand at room temperature and then invert carefully for 15 minutes to ensure homogeneity (rotation mixer).

Limitation of procedure

Discard a reference material if showing evidence of an insoluble sediment or microbiological contamination.

Attention

Each serum or blood donor unit used to manufacture this reference material was tested by state-accepted testing methods and found non-reactive for the following Antigens and Antibodies: HBsAG, anti-HIV-1, anti-HIV-2, anti-HBc and anti-HCV. No test method can offer complete assurance that products containing human source materials will be absent of these or other infectious agents. Therefore, handle them with the same precautions used with patient specimens.

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Medichem Diagnostica GmbH