Medidrug® Ethanol S-X X 1

 8°C 2°C	LOT AE03810  2027-03
Edition 16 / 2025-03	

Ethanol in Serum (human)

Werte-Ermittlung: ETOH 04/08

Ethanol in Serum (human)

Determination of values: ETOH 04/08

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Analysenmethode Method of Analysis Méthode d'analyse	Einheit Unit Unité	X 1
GC Gaschromatographie / Gas chromatography / Chromatographie gazeuse	g/L mmol/L	0,5 – 0,6 10,85 – 13,02

Anmerkungen

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In vitro diagnosticumEthanol
in Serum (human)**Description**

Medidrug® Ethanol S-X is a liquid reference material, prepared from 100% human serum in three significant levels.

Intended Use

Precision control of determinations of ethanol in serum.

Assignment of Values

Before using a new batch, Medidrug® Ethanol S-X first should be analyzed in a sufficient period of time. The assay values received apply only to the particular batch (Lot No. see label). Then the control ranges had to be calculated in accordance with relevant guidelines (see bibliography).

Transport and Storage

The reference material can be shipped at ambient temperature. Upon receipt the material has to be stored at +2 to +8°C.

Stability

Originally closed: +2 to +8°C, 36 months
But not longer than the expiration date specified in the certificate of analysis.
After working close the vials immediately and store them tightly capped at +2 to +8°C, do not let them standing open!
Ampoules are only for working at once!

Procedure

The reference material is ready to use and should be treated the same as patient specimens and run in accordance with the instructions accompanying the instrument.
Before sampling, allow the reference material to stand at room temperature and gently invert the vial to ensure homogeneity (5 min.).

Limitation of procedure

If there is evidence of microbial contamination or excessive turbidity in the used control, discard the vial.

Attention

Each serum or blood donor unit used to manufacture this control was tested by state-accepted testing methods and found non-reactive for following Antigens and Antibodies: HBsAG, anti-HIV-1, anti-HIV-2, anti-HBc and anti-HCV. No test method can offer complete assurance that products containing human source materials will be absent of these are other infectious agents. Therefore, handle them with the same precautions used with patient specimens.

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MEDICHEM Diagnostica GmbH