Edition 5 / 2024-11

**Alkoholkonsummarker
in Urin**

Werte-Ermittlung: ALCM RV 19-14

**Alcohol Consumption Markers
in Urine**

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6) Abweichung des Füllvolumens: max. ± 0,5%

Notes

1) The target value is calculated as the average of multiple analyses by independent laboratories accredited according to DIN EN ISO 17025 resp. DIN EN ISO 15189.
The analyses were performed by (UP)LC-MS-MS according to the guidelines of the German Medical Association for Quality Assurance (RiliBÄK) and the GTFCh (Society of Toxicological and Forensic Chemistry).

2) The estimated expanded uncertainty U of the target value is calculated with a coverage factor of k=2 corresponding to a level of confidence of about 95% as defined in the GUM.

3) Confidence interval:
Range in which the real value is located with a significance level of 95% (target value ± U, k=2).

4) The control ranges correspond to the accepted ranges used in the proficiency tests of the GTFCh. They are calculated from the target value and the double standard deviation according to Horwitz ($SD_{Horwitz}$): target value ± 2SD_{Horwitz} (last update 2020-01).
Regardless, the control ranges should be defined by the laboratory according to its own analytical requirements.

5) Matrix characteristics of this batch to assess the specimen validity according to SAMHSA guidelines:
pH: 7.3
Creatinine: 46.0 mg/dL
Density (20°C): 1.0070 ± 0.0001 g/cm³

6) Variation of the filling volume: max. ± 0.5%

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HORWITZ, W. (1982): Evaluation of analytical methods used for regulation of foods and drugs. Anal. Chem. 54 (1): 67A-76A

Mandatory Guidelines for Federal Workplace Drug Testing Programs, Federal Register January, 2017 (Vol. 82, No 13); SAMHSA Substance Abuse and Mental Health Services Administration

Guide to the Expression of Uncertainty in Measurement, ISO, 1995

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In-vitro-Diagnosticum

Alcohol Consumption Markers
in Urine

Description
Medidrug® ALCM U is a lyophilized reference material, prepared from human urine.

Intended use
The reference material is used for quality assurance of the determination of ethyl glucuronide and ethyl sulfate in urine.

Transport and Storage
The reference material can be shipped at ambient temperature. Upon receipt the material has to be stored at +2 to +8°C.

Stability
Originally closed: +2 to +8°C, 36 months
But not longer than the expiration date specified in the certificate of analysis.

After reconstitution: +2 to +8°C, 7 days

Preparation
Reconstitute each vial exactly with 2.5 ml of bi-distilled or HPLC-water. Replace the stopper, swirl gently and allow the reference material to stand for 30 minutes at room temperature, then gently invert the vial for 15 minutes (rotation mixer) to ensure homogeneity. After reconstitution the reference material should be treated the same as patient specimens and run in accordance with the instructions accompanying the instrument, kit or reagent being used. Before sampling, allow the reference material to sit at room temperature and gently invert the vial for 10 minutes (rotation mixer) to ensure homogeneity.

Precautions
No test method can guarantee complete assurance that products containing human source materials will be absent of these or other infectious agents. Therefore, handle them with the same precautions used with patient specimens.

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