Edition 5 / 2025-08

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3) Dichte (20°C): 0,9965 ± 0,0001 g/cm³

Notes

1) The verified value is calculated as the average of multiple GC analyses by laboratories accredited according to EN ISO 17025 resp. DIN EN ISO 15189.

2) The estimated expanded uncertainty U of the target value is calculated with a coverage factor of k=2 corresponding to a level of confidence of about 95% as defined in the GUM.

Confidence interval – statistic:
Range in which the real value is located with a significance level of 95% (target value ± U, k=2).

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Guide to the Expression of Uncertainty in Measurement, ISO, 1995

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Aqueous Standard

Description

Medidrug® BGS W high is an aqueous standard with a high concentration range of analytes.

Intended Use

For calibration or quality control of congener alcohols determination without dilution or diluted with blank serum.

Assignment of Values

The verified value is calculated as the average of multiple analyses by independent, accredited institutions of forensic medicine.

The stated values apply only to the particular batch (LOT see label).

Transport and Storage

The reference material can be shipped at ambient temperature. Upon receipt the material has to be stored at +20 to +25°C.

Stability

Originally closed and stored at +20 to +25°C the reference material is stable for 36 months, but not longer than the expiration date specified in the certificate of analysis.

Ampoules are only for working at once!

Procedure

Medidrug® BGS W high is ready to use and should be treated the same as patient specimens and run in accordance with the instructions accompanying the instrument, kit or reagent being used.

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MEDICHEM Diagnostica GmbH